

bs-18523R

[Primary Antibody]

ZFHX4 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Cow, Dog, Horse) Predicted MW.: 394 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 79776	SWISS: Q86UP3	
Target: ZFHX4		
Immunogen: KLH conjugated synthetic peptide derived from human ZFHX4: 2801-3000/3567.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: ZFHX4 (Zinc finger homeobox 4) is a putative transcription factor. It may have a role in neural and muscle differentiation.		