

D-DIMER Mouse mAb

Cat. No:	V1102
Name:	D-DIMER Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	E1D10
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted:Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	10mM PBS, pH7.4.
Applications:	ELISA
Antigen:	Native D-dimer protein derived from human plasma (bs-41172P)
Pairing Information:	

ELISA		Detection	
		V1101	V1102
Capture	V1101	-	-
	V1102	+++	-