

## D-Dimer Mouse mAb

Cat. No: V1102  
Name: D-Dimer Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: E1D10  
Host/Source: Mouse  
Subtype: IgG2b  
Reactivity: (predicted:&nbsp;Human)  
Format: Liquid  
Concentration:  $\geq 1\text{mg/ml}$   
Buffer: 10mM PBS, pH7.4.  
Applications: ELISA  
Antigen: Native D-dimer protein derived from human plasma ([bs-41172P](#))  
Pairing Information:

ELISA		Detection	
		V1101	V1102
Capture	V1101	-	-
	V1102	+++	-