

**bs-4649R**

**[ Primary Antibody ]**

## Parvovirus Capsid protein VP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Parvovirus Capsid protein VP1

**Immunogen:** KLH conjugated synthetic peptide derived from Hamster parvovirus H1 Capsid protein VP1: 501-600/734.

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Applications:** ELISA (1:5000-10000)

**Reactivity:** (predicted: Murine minute virus (strain MVM prototype), Hamster parvovirus H1)

**Predicted MW.:** 81 kDa