

**bs-10469R**

**[Primary Antibody]**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## Sulfadiazine Rabbit pAb

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> Sulfadiazine</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p>	<p><b>Isotype:</b> IgG</p>	<p><b>Applications:</b> ELISA (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Sulfadiazine)</p> <p><b>Predicted MW:</b> 0.25 kDa</p>
---	----------------------------	---

### — VALIDATION IMAGES —