

bs-7768R

[Primary Antibody]

RCC2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Chicken, Dog) Predicted MW.: 56 kDa Subcellular Location: Cytoplasm ,Nucleus
Clonality: Polyclonal		
GeneID: 55920	SWISS: Q9P258	
Target: RCC2		
Immunogen: KLH conjugated synthetic peptide derived from human RCC2: 441-522/522.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Required for completion of mitosis and cytokinesis. May function as a guanine nucleotide exchange factor for the small GTPase RAC1.		

— SELECTED CITATIONS —

- **[IF=2.634]** Zhang S et al.miR-196a Promotes Proliferation and Inhibits Apoptosis of Immature Porcine Sertoli Cells.(2019) DNA Cell Biol.Jan;38(1):41-48. WB ;Pig. 30427699