

**bs-8239R**

**[ Primary Antibody ]**

## FRMD7 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>IHC-P</b> (1:100-500) <b>IHC-F</b> (1:100-500) <b>IF</b> (1:100-500) <b>ELISA</b> (1:5000-10000)  <b>Reactivity:</b> (predicted: Human, Mouse, Rat, Pig, Cow, Dog, Horse)  <b>Predicted MW.:</b> 82 kDa  <b>Subcellular Location:</b> Secreted ,Extracellular matrix ,Cytoplasm ,Nucleus
<b>Clonality:</b> Polyclonal		
<b>GeneID:</b> 90167	<b>SWISS:</b> Q6ZUT3	
<b>Target:</b> FRMD7		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human FRMD7: 551-660/714.		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> Mutations in this gene are associated with X-linked congenital nystagmus. [provided by RefSeq, Dec 2008]		