

---

## ZNRF4 Rabbit pAb

Catalog Number: bs-9142R

Target Protein: ZNRF4

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:50-200), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse)

Predicted MW: 44 kDa

Entrez Gene: 148066

Swiss Prot: Q8WWF5

Source: KLH conjugated synthetic peptide derived from human ZNRF4: 131-230/429.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNRF4 contains 1 PA (protease associated) domain and 1 RING-type zinc finger.