bs-9142R

[Primary Antibody]

www.bioss.com.cn

sales@bioss.com.cn techsupport@bioss.com.cn

400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

ZNRF4 Rabbit pAb

GenelD: 148066 SWISS: Q8WWF5

Target: ZNRF4

Immunogen: KLH conjugated synthetic peptide derived from human ZNRF4:

131-230/429. < Extracellular >

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: ZNRF4 contains 1 PA (protease associated) domain and 1 RING-

type zinc finger.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit, Pig, Sheep, Cow, Chicken, Dog, Horse)

Predicted MW.: 44 kDa

Subcellular Location: Cell membrane