

---

## ZNF258 Rabbit pAb

Catalog Number: bs-7200R

Target Protein: ZNF258

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 148 kDa

Entrez Gene: 9204

Swiss Prot: O95789

Source: KLH conjugated synthetic peptide derived from human ZNF258: 101-200/1325.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF258 is part of a novel, putative, zinc binding motif (MYM) family which encodes proteins that maintain the repeats of the MYM motif.