

bs-7200R

[Primary Antibody]

ZNF258 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|---|--|--|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 9204</p> <p>Target: ZNF258</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ZNF258: 101-200/1325.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: ZNF258 is part of a novel, putative, zinc binding motif (MYM) family which encodes proteins that maintain the repeats of the MYM motif.</p> | <p>Isotype: IgG</p> <p>SWISS: O95789</p> | <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 148 kDa</p> <p>Subcellular Location: Nucleus</p> |
|---|--|--|