

---

## ZNF276 Rabbit pAb

Catalog Number: bs-7206R

Target Protein: ZNF276

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Zebrafish, Chicken, Dog, Horse)

Predicted MW: 67 kDa

Entrez Gene: 92822

Swiss Prot: Q8N554

Source: KLH conjugated synthetic peptide derived from human ZNF276: 401-500/614.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF276 may be involved in transcriptional regulation.