
ZNF282 Rabbit pAb

Catalog Number: bs-7207R

Target Protein: ZNF282

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Cow, Dog)

Predicted MW: 72 kDa

Entrez Gene: 8427

Swiss Prot: Q9UDV7

Source: KLH conjugated synthetic peptide derived from human ZNF282: 181-280/671.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF282 contains 1 KRAB domain and 5 C2H2-type zinc fingers. It belongs to the krueppel C2H2-type zinc-finger protein family and binds to the U5 repressive element (U5RE) of the human T cell leukemia virus type I long terminal repeat.