### bs-8231R

### [ Primary Antibody ]

# BIOSS ANTIBODIES www.bioss.com.cn

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn

## FAM104B Rabbit pAb

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 90268 SWISS: Q96BN8

Target: FAM104B

Immunogen: KLH conjugated synthetic peptide derived from human FAM104B:

21-115/115.

**Purification:** affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: With 181 million base pairs encoding around 1,000 genes,

chromosome 5 is about 6% of human genomic DNA. It is associated with Cockayne syndrome through the ERCC8 gene and familial adenomatous polyposis through the adenomatous polyposis coli (APC) tumor suppressor gene. Treacher Collins syndrome is also chromosome 5 associated and is caused by insertions or deletions within the TCOF1 gene. Deletion of the p arm of chromosome 5 leads to Cri du chat syndrome. Deletion of 5q or chromosome 5 altogether is common in therapy-related acute myelogenous leukemias and myelodysplastic syndrome. The FAM105B gene product has been provisionally designated FAM105B pending

further characterization.

**Applications: WB** (1:500-2000)

400-901-9800

Reactivity: Human

Predicted MW.: 13 kDa

Subcellular Location: Cytoplasm

#### VALIDATION IMAGES



Sample: PC-3 Cell (Human) Lysate at 40 ug Primary: Anti-FAM104B (bs-8231R) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 13 kD Observed band size: 13 kD