

bs-7888R

[Primary Antibody]

Cyclin L Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 57018</p> <p>Target: Cyclin L</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human Cyclin L/CCNL1: 351-450/526.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: CCNYL1 belongs to the cyclin family, Cyclin Y subfamily and contains 1 cyclin N terminal domain. The specific function of CCNYL1 is not yet known.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8N7R7</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Chimpanzee)</p> <p>Predicted MW.: 60 kDa</p> <p>Subcellular Location: Secreted ,Extracellular matrix ,Cell membrane ,Cytoplasm ,Nucleus</p>
--	--	---