

bs-5138R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

LZTFL1 Rabbit pAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) Reactivity: Rat (predicted: Human, Mouse, Rabbit, Cow, Chicken, Dog, Horse) Predicted MW.: 35 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 54585	SWISS: Q9NQ48	
Target: LZTFL1		
Immunogen: KLH conjugated synthetic peptide derived from human LZTFL1: 101-200/299.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: LZTFL1 (Leucine zipper transcription factor-like 1) contains a leucine zipper pattern and coiled-coil domains.		

— VALIDATION IMAGES —

Tissue/cell: rat spleen tissue; 4% Paraformaldehyde-fixed and paraffin-embedded; Antigen retrieval: citrate buffer (0.01M, pH 6.0), Boiling bathing for 15min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30min; Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min; Incubation: Anti-LZTFL1 Polyclonal Antibody, Unconjugated(bs-5138R) 1:200, overnight at 4°C, followed by conjugation to the secondary antibody(SP-0023) and DAB(C-0010) staining