

## T3 conjugated BSA

**Catalog Number :** bs-0339PB


**Quantity size:** 1mg

**Background:** T3(Triiodothyroxine)although only a minority of the total thyroxine, but the T3 is very important to check thyroid function. In the biochemical function of T4 than strong, and in the case of hyperthyroidism, T3 concentration in the blood are often marked increase (or even up to 5 times the normal).

CAS Number: 6893-02-3

Chemical structure:  $C_{15}H_{12}I_3NO_4$

MW: 650.97 g/mol



**Carrier protein:** BSA

**Storage:** Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

**Application Range:** Could use for immuring animal to prepare antibody or coating antigen to detect the antibody.

**Important Note:** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.