

Melamine conjugated KLH


Catalog Number : bs-2182PK

Quantity size: 1mg

Concentration: 1.25mg/ml, in 0.01M PBS(pH7.4)

Background: Melamine, also known as cyanuramide, is a synthetic chemical product that forms hard resins when combined with formaldehyde. It is used in a wide range of products such as cooking utensils, plates, industrial coatings, paper and paperboard, and flame retardant.

Chemical structure: $C_3N_6H_3$; MW: 126.12 g/mol



Carrier protein: KLH

Storage: Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

Application Range: Could use for immunizing animal to prepare antibody or coating antigen to detect the antibody.

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.