

Mouse IgG2b, Isotype Control

Cat. Number:	bs-30102P
Purity:	≥95%, Purified by Protein G
Form:	Liquid
Concentration:	1 mg/ml.
Buffer:	0.01M PBS(pH7.0)
Storage:	Store at -20℃ for one year. Avoid repeated freeze/thaw cycles.
Application:	Isotype Control 0.125ug/1X10 ⁶ cells
Usage:	The reagent is intended as isotype control for flow cytometry analysis to establish the amount of non-specific antibody binding. For your particular experiment, use the same concentration of this isotype control antibody as the recommended working concentration of the antigen-specific antibody. If under particular experimental conditions the background signal of the isotype control is too high (usually when working concentrations of used antibodies are above 10 µg per ml of incubation mixture), change the conditions of your experiment to reduce the background.