

Recombinant Dengue protein (contains 13 immuno-reactive epitpes)

Cat. Number:	bs-41421P
Product name:	Recombinant Dengue protein (contains 13 immuno-reactive epitpes)
Expression system:	E.coli
Tags:	Tag free
Detected MW:	16 kDa
Format:	Liquid
Purity:	>95% as determined by SDS-PAGE
Concentration:	>0.5 mg/ml
Buffer:	PBS (pH=7.4)
Storage:	Stored at -70℃ or -20℃. Avoid repeated freeze/thaw cycles.
Important Note:	This product as supplied is intended for research use only, not for use in human,therapeutic or diagnostic applications.
Description:	This protein is a dengue protein for dengue antibody rapid test. This protein contains 13 immuno-reactive epitpes derived from dengue genome,it reacts to most four serotype dengue antibodies.The rapid test prepared by this protein has over 90% sensitivity and specifity over 95% to detect both dengue IgM and IgG.