


Product Details

Product name:	Anti-human IL-6Ra / CD126 (vobarilizumab Biosimilar)	SKU:	BIO0044SM
Target Name:	IL-6Ra / CD126	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P08887	Concentration:	Lyophilized
Clone#:	vobarilizumab	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	145 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-IL-6Ra / CD126 (vobarilizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 96.9%.

Purity:SEC-HPLC




The purity of Anti-IL-6Ra / CD126(vobarilizumab) is 99.34%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human IL 6R His at 2 ug/mL can bind Anti-IL-6Ra / CD126 (vobarilizumab), EC50=0.0158 ug/mL.

Bioactivity:FACS



Human IL6R CHO cells were stained with Anti-IL-6Ra / CD126 (vobarilizumab) and negative control protein respectively, washed and then followed by PE and analyzed with FACS, EC50=1.0090 nM.