


Product Details

Product name:	Anti-human Siglec-15 / CD33L3 (NC318 Biosimilar)	SKU:	BIO0073SM
Target Name:	Siglec-15 / CD33L3	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q6ZMC9	Concentration:	Lyophilized
Clone#:	NC318	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	146.98 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE



Anti-Siglec-15 / CD33L3 (NC318) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-Siglec-15 / CD33L3 (NC318) is 100%, determined by SEC-HPLC.