


Product Details

Product name:	Anti-human ANGPTL3 (evinacumab Biosimilar)	SKU:	BIO0092SM
Target Name:	ANGPTL3	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q9Y5C1	Concentration:	Lyophilized
Clone#:	evinacumab	Isotype:	Human IgG4SP
Reactivity:	Human, Mouse, Rat, Cynomolgus, Rhesus	Calculated M.W.:	150 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE



Anti-ANGPTL3 (evinacumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-ANGPTL3 (evinacumab) is 99.48%, determined by SEC-HPLC.