


Product Details

Product name:	Anti-human CXCR5 / CD185 (SAR113244 Biosimilar)	SKU:	BIO0146SM
Target Name:	CXCR5 / CD185 / BLR1 / MDR15	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P32302	Concentration:	Lyophilized
Clone#:	SAR113244	Isotype:	Human IgG4SP
Reactivity:	Human	Calculated M.W.:	144.1 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE



Anti-CXCR5 / CD185 (SAR113244) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 98.9%.

Purity:SEC-HPLC



The purity of Anti-CXCR5 / CD185(SAR113244) is 98.97%, determined by SEC-HPLC.