


Product Details

Product name:	Anti-human LILRB2 / ILT4 / CD85d (JTX-8064 Biosimilar)	SKU:	BIO0221SM
Target Name:	LILRB2 / ILT4 / CD85d	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q8N423	Concentration:	Lyophilized
Clone#:	JTX-8064	Isotype:	Human IgG4
Reactivity:	Human	Calculated M.W.:	146.48 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-LILRB2 / ILT4 / CD85d (JTX-8064) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Purity:SEC-HPLC



The purity of Anti-LILRB2 / ILT4 / CD85d(JTX-8064) is 98.96%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human LILRB2 His at 2 µg/mL can bind Anti-LILRB2 / ILT4 / CD85d (JTX-8064), EC50=0.01229 µg/mL.