


Product Details

Product name:	Anti-human PCSK9 (alirocumab Biosimilar)	SKU:	BIO0245SM
Target Name:	PCSK9	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q8NBP7	Concentration:	Lyophilized
Clone#:	alirocumab	Isotype:	Human IgG1
Reactivity:	Human, Mouse, Cynomolgus	Calculated M.W.:	146.24 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-PCSK9 (alirocumab) on SDS-PAGE underreducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-PCSK9 (alirocumab) is 96.96%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human PCSK9 chis at 2 ug/mL can bind Anti-PCSK9 (alirocumab), EC50=0.04823 ug/mL.