


## Product Details

<b>Product name:</b>	Anti-human TNFSF11 / RANKL / CD254 (denosumab Biosimilar)	<b>SKU:</b>	BIO0255SM
<b>Target Name:</b>	TNFSF11 / RANKL / CD254	<b>Size:</b>	1 mg/5 mg/20 mg
<b>Target Uniprot:</b>	O14788	<b>Concentration:</b>	Lyophilized
<b>Clone#:</b>	denosumab	<b>Isotype:</b>	Human IgG1
<b>Reactivity:</b>	Human, Mouse, Rat	<b>Calculated M.W.:</b>	145.16 kDa
<b>Application:</b>	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	<b>Endotoxin:</b>	<0.001 EU/ug
<b>Formulation:</b>	0.1M Pro, 20mM Arg, pH5.0	<b>Conjugation:</b>	None
<b>Storage:</b>	For long term storage, the produce should be stored at -20° C or lower.	<b>Expression System:</b>	CHO
<b>Reconstitution:</b>	Dissolve with sterile ddH2O	<b>Purification:</b>	Protein A


## Data

### Purity:SDS-PAGE



Anti-TNFSF11 / RANKL / CD254 (denosumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

### Purity:SEC-HPLC



The purity of Anti-TNFSF11 / RANKL / CD254(denosumab) is 97.25%, determined by SEC-HPLC.