


## Product Details

Product name:	Anti-human Tau (gosuranemab Biosimilar)	SKU:	BIO0266SM
Target Name:	Tau	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P10636	Concentration:	Lyophilized
Clone#:	gosuranemab	Isotype:	Human IgG4SP
Reactivity:	Human	Calculated M.W.:	144.94 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


## Data

### Purity:SDS-PAGE




Anti-Tau (gosuranemab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

### Purity:SEC-HPLC



The purity of Anti-Tau (gosuranemab) is 97.31%, determined by SEC-HPLC.

### ELISA



Immobilized human Tau 441, His Tag at 2 g/mL can bind Anti-Tau (gosuranemab) EC<sub>50</sub>=0.004731 ug/mL.