


Product Details

| | | | |
|------------------------|---|---------------------------|-----------------|
| Product name: | Anti-human GFRAL (NGM120 Biosimilar) | SKU: | BIO0308SM |
| Target Name: | GFRAL | Size: | 1 mg/5 mg/20 mg |
| Target Uniprot: | Q6UXV0 | Concentration: | Lyophilized |
| Clone#: | NGM120 | Isotype: | Human IgG1 |
| Reactivity: | Human | Calculated M.W.: | 145.7 kDa |
| Application: | ELISA, Bioactivity: FACS, Functional assay, Research in vivo | Endotoxin: | <0.001 EU/ug |
| Formulation: | 0.1M Pro, 20mM Arg, pH5.0 | Conjugation: | None |
| Storage: | For long term storage, the produce should be stored at -20° C or lower. | Expression System: | CHO |
| Reconstitution: | Dissolve with sterile ddH2O | Purification: | Protein A |


Data

Purity:SDS-PAGE



Anti-GFRAL (NGM120) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-GFRAL (NGM120) is 100%,determined by SEC-HPLC.