


## Product Details

Product name:	Anti-human CD28 (lulizumab Biosimilar)	SKU:	BIO0357SM
Target Name:	CD28	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Phage display	Concentration:	Lyophilized
Clone#:	lulizumab	Isotype:	Human V-kappa
Reactivity:	Human,Cynomolgus [US20110097339A1]	Calculated M.W.:	12.01 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


## Data

### Purity:SDS-PAGE



Anti-CD28 (lulizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

### SEC-HPLC



The purity of Anti-CD28 (lulizumab),determined by HPLC, is 100%