


Product Details

Product name:	Anti-human MRC2 / CD280 (Copenhagen Rigshospitalet patent Anti-human uPARAP Biosimilar)	SKU:	BIO0363SM
Target Name:	MRC2 / CD280	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q9UBG0	Concentration:	Lyophilized
Clone#:	Copenhagen Rigshospitalet patent Anti-human uPARAP	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	145.28 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-MRC2 / CD280 (Copenhagen Rigshospitalet patent anti-uPARAP) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Purity:SEC-HPLC



The purity of Anti-MRC2 / CD280 (Copenhagen Rigshospitalet patent anti-uPARAP) is more than 95%, determined by SEC-HPLC.

ELISA



Immobilized Human MRC2 / CD280, His Tag can bind Anti-MRC2 / CD280 (Copenhagen Rigshospitalet patent anti-uPARAP). The EC50 was approximately 0.01620 µg/mL.