


Product Details

Product name:	Anti-human TNFRSF5 / CD40 (selicrelumab Biosimilar)	SKU:	BIO0404SM
Target Name:	TNFRSF5 / CD40	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P25942	Concentration:	Lyophilized
Clone#:	selicrelumab	Isotype:	Human IgG2SA
Reactivity:	Human	Calculated M.W.:	146.2 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-TNFRSF5 / CD40 (selicrelumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-TNFRSF5 / CD40 (selicrelumab) is more than 95%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human CD40 His at 2 ug/mL can bind Anti-TNFRSF5 / CD40 (selicrelumab), EC50=0.009782 ug/mL.