


Product Details

Product name:	Anti-human B7-H2 / ICOSL / CD275 (prezalumab Biosimilar)	SKU:	BIO0438SM
Target Name:	B7-H2 / ICOSL / CD275	Size:	1 mg/5 mg/20 mg
Target Uniprot:	O75144	Concentration:	Lyophilized
Clone#:	prezalumab	Isotype:	Human IgG2SA
Reactivity:	Human	Calculated M.W.:	145 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-B7-H2 / ICOSL / CD275 (prezalumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-B7-H2 / ICOSL / CD275 (prezalumab) is 98.05%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human B7 H2 His at 2 ug/mL can bind Anti-B7-H2 / ICOSL / CD275 (prezalumab), EC50=0.007111 ug/mL.