


Product Details

Product name:	Anti-human CEACAM1 / CD66a (Chaim Sheba Med. Cntr. patent Anti-human CEACAM1 Biosimilar)	SKU:	BIO0506SM
Target Name:	CEACAM1 / CD66a	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P13688	Concentration:	Lyophilized
Clone#:	Chaim Sheba Med. Cntr. patent Anti-human CEACAM1	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	144.98 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE



Anti-CEACAM1 / CD66a (Chaim Sheba Med.Cntr. patent anti-CEACAM1) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-CEACAM1 / CD66a (ChaimSheba Med. Cntr. patent anti-CEACAM1) is more than 95%, determined by SEC-HPLC.