


Product Details

Product name:	Anti-human Amyloid Beta (Georgetown U. patent Anti-human Amyloid Beta Biosimilar)	SKU:	BIO0561SM
Target Name:	Amyloid Beta	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P05067	Concentration:	Lyophilized
Clone#:	Georgetown U. patent Anti-human Amyloid Beta	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	150 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-Amyloid Beta (Georgetown U. patent anti-Amyloid Beta) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Purity:SEC-HPLC



The purity of Anti-Amyloid Beta (Georgetown U. patent anti-Amyloid Beta) is more than 95%, determined by SEC-HPLC.

ELISA



Immobilized Human Amyloid Beta1-42, No Tag at 2 ug/mL can bind Anti-Amyloid Beta (Georgetown U. patent anti-Amyloid Beta), EC50=0.005037 ug/mL.