


Product Details

Product name:	Anti-human Integrin αV / ITGAV / CD51 (abrituzumab Biosimilar)	SKU:	BIO0688SM
Target Name:	Integrin αV / ITGAV / CD51	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P06756	Concentration:	Lyophilized
Clone#:	abrituzumab	Isotype:	Human IgG2SA
Reactivity:	Human	Calculated M.W.:	144.58 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE



Anti-Integrin αV / ITGAV / CD51 (abrituzumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-Integrin αV / ITGAV / CD51 (abrituzumab) is 98.44%, determined by SEC-HPLC.