


## Product Details

Product name:	IgG1/Kappa Isotype Control	SKU:	BIO0701SM
Target Name:	Spike RBD	Size:	1 mg/5 mg/20 mg
Target Uniprot:	None	Concentration:	Lyophilized
Clone#:	Monoclonal	Isotype:	Human IgG1
Reactivity:	Virus	Calculated M.W.:	150 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


## Data

### FACS




This IgG1/kappa isotype antibody can be used in binding Assay with Pan-T cell Line as a antibody negative control

### FACS




This IgG1/kappa isotype antibody can be used in binding Assay with CHOS as a antibody negative control

### FACS



This IgG1/kappa isotype antibody can be used in binding Assay with HEK293 as a antibody negative control

### FACS



This IgG1/kappa isotype antibody can be used in binding Assay with NCI-H226 as a antibody negative control