


Product Details

Product name:	Anti-human Nectin-4 (enfortumab vedotin-ejfv Biosimilar)	SKU:	BIO0863SM
Target Name:	Nectin-4	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q96NY8	Concentration:	Lyophilized
Clone#:	enfortumab vedotin-ejfv	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	143.9 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	MMAE(ADC)
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-Nectin-4 (enfortumab vedotin-ejfv) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-Nectin-4 (enfortumab vedotin-ejfv) is 100%, determined by SEC-HPLC.

Bioactivity:ELISA



Immobilized human NECTIN 4 His at 2 µg/mL can bind Anti-Nectin-4 (enfortumab vedotin-ejfv), EC₅₀=0.003846 µg/mL.