


Product Details

Product name:	Anti-human ALCAM/CD166 (praluzatamab ravtansine Biosimilar)	SKU:	BIO0883SM
Target Name:	ALCAM / CD166	Size:	1 mg/5 mg/20 mg
Target Uniprot:	Q13740	Concentration:	Lyophilized
Clone#:	praluzatamab ravtansine	Isotype:	Human IgG1
Reactivity:	Human,Cynomolgus [US10745481]	Calculated M.W.:	145 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	DM4(ADC)
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity: SDS-PAGE



Anti-ALCAM/CD166 (praluzatamab ravtansine) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Purity: SEC-HPLC



The purity of Anti-ALCAM/CD166 (praluzatamab ravtansine) is more than 95%, determined by SEC-HPLC.