

Product Details

Product name:	Anti-human CEACAM5 / CEA / CD66e (Tusamitamab ravtansine Biosimilar)	SKU:	BIO0885SM
Target Name:	CEACAM5 / CEA / CD66e	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P06731	Concentration:	Lyophilized
Clone#:	Tusamitamab ravtansine	Isotype:	Human IgG1
Reactivity:	Human,Cynomolgus [US9617345]	Calculated M.W.:	144.52 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	DM4(ADC)
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A

Data