


Product Details

Product name:	Anti-human TM4SF1 (Beth Israel Patent Anti-human TM4SF1 Biosimilar)	SKU:	BIO0917SM
Target Name:	TM4SF1 / T4S1 / M3S1 / TAAL6	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P30408	Concentration:	Lyophilized
Clone#:	Beth Israel Patent Anti-human TM4SF1	Isotype:	Human IgG1
Reactivity:	Human	Calculated M.W.:	150 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.1M Pro, 20mM Arg, pH5.0	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity: SDS-PAGE



Anti-TM4SF1 (Beth Israel Patent Anti-TM4SF1) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.

Purity: SEC-HPLC



The purity of Anti-TM4SF1 (Beth Israel Patent Anti-TM4SF1) is more than 95%, determined by SEC-HPLC.