


Product Details

Product name:	Anti-human HGFR / c-Met (amivantamab Biosimilar)	SKU:	BIO0014SM
Target Name:	HGFR / c-Met	Size:	1 mg/5 mg/20 mg
Target Uniprot:	P00533 & P08581	Concentration:	Lyophilized
Clone#:	amivantamab	Isotype:	Human IgG1-kappa
Reactivity:	Human	Calculated M.W.:	145.88 kDa
Application:	ELISA, Bioactivity: FACS, Functional assay, Research in vivo	Endotoxin:	<0.001 EU/ug
Formulation:	0.01M PBS, pH7.4.	Conjugation:	None
Storage:	For long term storage, the produce should be stored at -20° C or lower.	Expression System:	CHO
Reconstitution:	Dissolve with sterile ddH2O	Purification:	Protein A


Data

Purity:SDS-PAGE




Anti-HGFR / c-Met (amivantamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Purity:SEC-HPLC



The purity of Anti-HGFR / c-Met (amivantamab) is more than 95%, determined by SEC-HPLC.

Bioactivity:ELISA



Detects EGFR/ERBB1/HER1 in indirect ELISAs.