



Human IgE, unconjugated

Cat. Number: bs-2226P

Background: Human IgE full length protein was purified by Protein L chromatography. As the IgE comes from a monoclonal cell line, there is no contamination of antibodies of other isotypes.

The cell line supernatant is non-reactive for HIV1/HIV2 antibodies, HBsAg, Hepatitis C antibodies, HIV1 and Hepatitis C by NAT as per FDA methods.

Application:

- ELISA: Use at an assay dependent dilution.
- Not yet tested in other applications.
- Optimal working dilutions must be determined by the end user.

Storage: Shipped at 4°C, Store below -20°C (Avoid repeated freeze/thaw cycles).

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.